




AURIGENE

PHARMACEUTICAL SERVICES



Case Study

How the re-designed synthesis of a complex carbohydrate can reduce cycle time.



Project challenge:

Complex carbohydrate chemistry involving a linear sequence of 10 chemical transformations, unstable intermediates and column chromatographic purifications

Solution design:

Process research and optimization was performed to develop a robust and scalable process which was implemented on commercial scale.

- Telescoping of reactions reduced the number of isolable intermediates thereby reducing the cycle time by 50% >
- Risk of handling highly unstable intermediates was mitigated.
- Use of Stereoselective glycosylation to maximize the desired stereoisomer
- Efficient crystallization procedure developed by identifying the onset conditions thereby avoiding column chromatography purification



Conclusion:

Development and manufacturing of carbohydrates bearing labile functional groups is exemplified with this example. Several molecules where protection-deprotection strategy and highly stereoselective transformations are involved have been worked on by our team.

Thank You



For more information please visit
<https://www.aurigeneservices.com/>



To place an inquiry please visit
<https://www.aurigeneservices.com/form/contact>



Mail us at
contactapsl@aurigeneservices.com