



# AURIGENE

PHARMACEUTICAL SERVICES



# Cell Engineering

Recombinant cell line development services for wide protein classes

Tool Reagents | Biologic Production


 /AurigenePharma

 /aurigenepharmaceuticalservices

 /AurigenePharmaceuticalServices

 [www.aurigeneservices.com](http://www.aurigeneservices.com)

 [contactapsl@aurigeneservices.com](mailto:contactapsl@aurigeneservices.com)


 North America: +1 617 821 0595

 European Union: +32 471 808710

# Why Aurigene Pharmaceutical Services?



Recombinant cell lines for wide protein classes viz. RTK, GPCR, transporters, ion channels, dual or co-receptors, immune modulation (ligand/receptors) and cell adhesion molecules



Complex target recombinant cell lines



Cell line with specific receptor numbers



Double or multiple gene KD / KO cell lines



Immune cells (T-cells / B-cells) reporter cell lines



In-house lentiviral platform for producing difficult to express reagent proteins

## Recombinant cell line types:

Membrane bound receptor cell lines | Knock-Down (KD) cell lines | Knock-Out (KO) cell lines | Inducible cell lines | Reporter cell lines | Protein production cell lines (reagent & therapeutic)

## Technologies:

Standard (Lipid / nanoparticle) | Electroporation | Transposon | Lentiviral / Retroviral | Tet or Cumate induction | CRISPR/Cas9 | TALEN | ZFN | siRNA / shRNA | Fluorescent protein and Luciferase

## Screening:

Recombinant cell lines: Flow cytometry | Western blot | qRT-PCR | Sequencing  
Protein production cell lines: ELISA | ClonePix | PA-HPLC | Octet

## Quality Checks:

### Recombinant cell lines:

Thaw test | Mycoplasma | Receptor density estimation | Clone stability analysis

### Protein production cell lines:

Clone stability analysis (30 generation) | Mycoplasma | Endotoxin | Bacterial contamination | SE-HPLC | LC-MS analysis Glycoform and charge variant analysis